

# DUAL-MODE Infrared Thermometer

## JPD-FR301



**JPD-FR301** Infrared Thermometer measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can quickly get measurement results after positioning properly the temperature probe in the ear canal or forehead.



### Magnetic Probe Cover

Automatically switched between forehead & ear mode

#### Safe & Hygienic



#### Multiple use




#### 20 memories



#### Fever warning





## SPECIFICATION

Display	HD LED Display
<b>Measuring range</b>	
Forehead mode	22-43°C / 71.6-109.4°F
Ear mode	34-43°C / 93.2-109.4°F
Object mode	0-100°C / 32-212°F
<b>Accuracy(Laboratory)</b>	
Forehead mode:	± 0.2 °C (36-39 °C); ± 0.3 °C (22-36 °C / 39-43 °C )
Ear mode:	± 0.2 °C (36-39 °C); ± 0.3 °C (34-36 °C / 39-43 °C )
Object mode:	± 1.0°C / ± 2.0°F
Battery	AAAx2
Memory storage	20 sets
Temperature units	°C/°F adjustable
Automatic shutdown function	10s ± 1s
Display resolution	0.1°C/°F

## STRUCTURE

- |                    |                             |
|--------------------|-----------------------------|
| (1) Probe cover    | (5) Battery cover           |
| (2) Probe          | (6) Adult/Child/Object mode |
| (3) LED display    | (7) Memory/Mute             |
| (4) Power/ Measure | (8) °C/°F switch            |



## PACKING SPECIFICATION

CTN SIZE(CM)	QTY/CTN	GW(KG)	FCL/20
54x29x40	100	20	400

### Accessories

AAA Batteries x 2 | Storage bag | User Manual | Quick guide x 1

